


Electrode to Connector Map		
Assembly type	Any 'C' style 32-channel probe on C32 package (standard).	
Version	V6. July 2007	
Connector	Two Omnetics nano NPD series (A7695-001)	
Mating Connector	Two Omnetics nano NPD series (A7519-001)	
Compatible Headstages	TDT	Two RA16CH, LP16CH
	Multi-Channel Systems	Two MPA8I
	Plexon	Two HST/16o25-18P, HST/16V-G20
	Neuralynx	Two HS-16
	Alpha Omega	Two 16 mini, TeleSpike
Compatible Adaptors	TDT	AC-CH
	Reference site	Present on some designs. Please check the catalog CAD drawings.
Note	Electrode sites face down (as show in the picture).	

CONNECTOR PIN OUT



Orientation: probe points to top as shown



ELECTRODE SITE MAP

